



CHECKLIST

(See page 3 for visual identification of checklist parts)

<u>Part:</u>	<u>Model number:</u>
__ Autoclave	
__ Door seal	200-A00003
__ Tray	100-A00104
__ Tray support	100-A00105
__ Manual	310-A00001
__ Power cord	400-P00004

Notes & Service Record

For your own future reference and to expedite answers to your questions, please complete the following information:

Date Purchased: _____

Serial Number: _____

Model Number: _____

Purchased From: _____

Your Notes:

Revolutionary Science recommends the weekly use of G. Bacillus

Stearotherophilus biological indicator ampules. Place indicator in center of full load and test weekly.

The use of FDA cleared class 5 integrators should be used with each load. Follow the manufacturer's recommendations for class 5 integrator use.